



### INFUSION BOMB ENTERAL DAIHWA DI-4000

MEDIFUSION can be used for continuous flow rate infusion of fluid drug and transfusion in ICU or CCU as well as infusion of anesthetic, anticancer agent or nutrient. With any standard IV set without special trying & efforts, the operator can set the pump for the application.

#### EXPLANATION

- Dural ability of docking station
- Stand-alone system
- Drug library
- IV set library max 300
- Anti-bolus
- Language configuration
- Secondary mode
- Easy user interface
- Feeding for enteral nutrition

#### STANDARD

1. Flow rate : 0.1~1200 ml/h
2. Infused Volume : 0.1~9999 ml
3. Bolus rate : 1~1200 ml/h
4. K.V.O : 0.1~9 ml/h
5. Occlusion Level : 9 Level (60-850 mmHg)
6. Battery : Operating 6 hour, Recharging 5 hour
7. Weight : 1.8 kg