



INFUSION PUMP 2 CHANNELS DAIHWA DI-2200

2-Channel Pump General Ward, ICU and CCU

EXPLANATION

1. Practical use in ICU
2. Two (2) separate CPU (Dual Channel)
3. Save Space & (Economic)
4. Enable to make operator's own IV set list in the software after calibration
5. Compact Design, lightweight, easy portable & safety proved
6. Improved Pumping unit – Accuracy improved

FEATURES

1. Compact Design
2. Small size – 3.0 kg
3. Compatible with various kinds of IV set
4. One touch Bolus, Pause function
5. Dosage – auto calculation (Input data: body weight, drug volume, solution volume, dose rate)
6. Titration – Enable to change the flowrate during infusion
7. Keylock for safety
8. Event history up to 2000
9. Firmware upgrade via internet (system or software upgrade)
10. Quality Assurance & Certification – ISO13485, CE

STANDARD

1. Flow Rate : 0.1 ~ 1200 ml/h
2. Accuracy : 5% with assigned IV set



- 3. Total Volume : 0 to 9,999ml or no limit
- 4. Bolus/Purge : 700ml/h(default)/ 1~1200ml/h
- 5. Battery : Charging – 10 hr , Operating – 4hr @ 100ml/h
- 6. Other Function
Titration, Date and Time, History Monitoring (up 1000)
Firmware Upgrade Program